



Medizintechnik Sattler GmbH

Innovations for the surgery

Systemdescription:

The „Sattler“ Königsee retractor system is a universally designed surgical retractor system. The system is modular in construction and can be adapted and expanded any time as desired. The self-retaining modular construction makes it efficient and optimally cost-effective to use. This can reduce operation costs so it is an investment that will soon pay for itself.

Areas of use:

- Abdominal
- Gynecology
- Urology
- Paediatric surgery
- Vascular surgery

Product features:

100% Made in Germany

broad range of uses

stable retaining system



SATTLER RETRACTOR SYSTEM

RETRACTOR-KIT/ ABDOMINAL VERSION



Pos.	Item number	Rilsan	articulation	description	Classic-Line	Premium-Line
1	01100-52			sterile field post, one clamp 20 / 14 mm	1	1
2	01115-14			support arm asymmetric, two clamps 14/14 mm	1	1
3	01120-14			frame arm (couple) for retractor blades	1	1
4	01130-05			excenter snap clamp 14 / 8 mm	5	7
5	0103	x		blade Fritsch 60 x 50 mm		2
6	0104	x		blade Fritsch 80 x 60 mm	2	2
7	0105	x		blade Fritsch 80 x 80 mm	2	2
8	0107	x		blade Fritsch 90 x 105 mm		2
9	0205	x	x	blade Langenbeck 30 x 190 mm		1
10	0210	x	x	blade Langenbeck 50 x 170 mm	1	1
11	0211	x	x	blade Langenbeck 50 x 220 mm	1	1
12	0213	x	x	blade Langenbeck 70 x 190 mm		1
13	0410	x	x	blade Deaver 50 x 170 mm		1
14	0413	x	x	blade Deaver 70 x 190 mm	1	1
15	0501	x		blade Deaver slotted 75 x 175 mm	1	1
16	0702	x		blade Harrington 70 x 150 mm	1	1
17	1101		x	blade finger fence 150 x 90 mm,flexible	1	1
18	1302	x	x	blade "Body wall" 80 x 160 mm	1	1
19	03FE000273			cold light cable Ø 4,5/4,8x2650 mm elastically flexible, with bladeholder		1
20	01130-12			universal retaining system for MIC- instruments		1